**UK text:**

**AIR 2:**

Break free from the endless cycle of buying razors, wax strips, and inconvenient trips to the clinic. Ulike IPL laser hair removal saves you time and budget for more important things in life. Unique Ice-cooling technology offers a cool touch while the device delivers efficient energy to target hair follicles to break the growth cycle and achieve a long-lasting result.

Comfort and effectiveness, you can have both. Keep your skin flawless and say goodbye to the last-minute hair removal hassle. You can get visible salon-quality results in 4 weeks\* (\*refers to a decrease in hair count within 4 weeks, sourced from third-party laboratory testing reports. Individual results may vary).

Our laser hair removal device is wired for 5 intensity levels available, from upper lips, underarms to thighs, you can always find the most efficient level for your needs, one device is suitable for whole body hair removal. With FDA clearance, expert backing, and over a decade of specialization in IPL skin safety, the Ulike Air 2 boasts a built-in UV filter along with our patented Sapphire Ice-Cooling Technology, ensuring gentle yet highly effective, nearly painless hair removal.

**Suitable for**

Adults

(Note: Ulike Air 2 is a wired device and is not suitable for tattooed and overly-tanned skin, wounded areas, and during pregnancy.)

**How to use**

We strongly suggest you complete a test patch with the lowest setting mode, note your skins reaction and proceed accordingly increasing the level of intensity as needed.

**Choose the Treatment Method**

- Auto-Glide Method - Press and hold for 2 seconds to activate, glide the light window on skin, the light will automatically flashes Great for large areas such as arms, legs
- Spot Method - short press the treatment button, one flash at a time, suitable for small areas such as fingers and upper lips

**Step 1 - Prepare:**

- Remove all visible hair by shaving
- Cleanse the treatment area and pat dry
- DON'T use any gels or creams with this device!

**Step 2 - Connect The Power:**

- Find a comfortable place
- Connect adapter to the handset
- Plug in the power

**Step 3 Goggle Up:**

- Wear the professional IPL Safety Glasses comes in the box

**Step 4 - Turn On The Device:**

- Short press to turn on the handset
- Once the handset is on, short press to change the intensity levels

**Step 5 - Full Contact:**

- Place the handset in a vertical position against skin.
- Make sure the skin sensor indicator is light up

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